

Adjustable Hand Crimper from Kebby Industries, Inc.

***Adjustable
Hand
Crimpers***



- Adjustable block corrects for differences in liner thicknesses

The Adjustable Crimper has an adjustable pressure block in the crimper head that allows crimping of liners with a thickness range of 0.25 to 3.0mm. The unit looks like the standard hand crimper but it is supplied with an allen wrench and special instructions for making adjustments.

Adjustable Crimper Instructions

The crimper has an adjustment feature which will allow it to be used with a variety of septa, discs, stoppers or lined caps. The user may adjust the stroke or movement for each application as needed.

Kebby Industries, Inc.
Rockford, IL
USA
Ph. 815-963-1466
sales@kebyindustries.com