

***Kebby Industries, Inc.***

***Micro-XL 4300  
Bench Crimper***

***Compact Series, with  
stainless-steel base***

***P/N A10021***



**Uses interchangeable crimping and decapping heads.  
(Sold separately)**

**Weight- 16 lbs**

**Unit Height- 14"**

**Base Dimensions- 8" x 10"**

**Bottle and cap assembly clearance height- 5"**

Kebby Industries, Inc.  
4075 Kilburn Ave.  
Rockford, IL 61101  
Ph. 815-963-1466  
Email [sales@kebyindustries.com](mailto:sales@kebyindustries.com)